





CLEANING OF MICROSCOPE OPTICS



GLYC (ne = 1.460

## www.leica-microsystems.de

Order no.: English 11914709 - Copyright © by Leica Microsystems CMS GmbH, Wetzlar, Germany, 2018. Subject to modifications. LEICA and the Leica Logo are registered trademarks of Leica Microsystems IR GmbH.

# 1 Considerations before cleaning

## Which image quality do you expect?

Optimize all microscope settings:

Aperture, Köhler, correction ring, etc.

## Check your specimen

I Use standard with proven quality

II Use cleaned specimen

## Check slide and cover glass

#### **Check immersion medium**

I Same type as embedding medium of specimen

II Is immersion medium from the same lot?

III Only use recommended types

IV No Anisol

# 2 Where to clean and where not



#### Where to clean

- 1 Front lens of objective
- 2 Cover glass of camera sensor
- 3 Cover slip and specimen slide
- 4 C-mount
- 5 Condenser lens
- 6 Misc. glass surfaces

#### Where not to clean

- I Filter cubes
- II All inside optics of a microscope

# 3 How to locate dirt

## Turn camera and inspect image

I If dirt is turning in the image  $\rightarrow$  not located in camera

II If dirt stays → camera

## Turn objective, condenser, C-mount

I If dirt is turning with image → objective, condenser, etc.

## Control the relevant elements macroscopically

I Magnifying glass or oculars upside down

#### Front lens

I Turn it

II Look from the back against a bright background

# 4 Different kinds of dirt

# Loosely or not permanently attached dirt

I Glass (broken slides, cover slips, etc.)

II Skin, dander

III Pollen, etc.

#### **Attached dirt**

I Water soluble

II Solvent soluble

III In practice often a mixture

# 5 How to clean

## Aim of cleaning

I Complete elimination of dust and dirt for perfect image quality

II No remaining residues on the optical parts

III No damage of optical parts

# **General procedure**

I Locate dirt

II Inspect dirt

III Remove dust

IV Remove water soluble dirt

V Remove solvent soluble dirt

VI Inspect result

VII Repeat from III to VI if necessary

### Remove dust

I If possible do it manually

II Only use specified compressed air

**a** Contamination with oil may cause difficulties to remove residues

III Perfect tool: bellows

## Remove attached dirt

I Never use rubbing materials, papers, microfibre cloths

II Never clean dry

III Always start to clean water soluble dirt

a If you can see grease start with solvents

IV All solvents for cleaning should be absolutely clean

a PA solvents (ultra pure)

**b** Destilled or demineralized water

**c** No ammonium containing glass cleaner like Sidolin, Sparkle, etc.

V All tools should be absolutely clean

a Use certified wound cotton sticks and dental cotton pads to avoid lint (no Q-tips or similar!)

#### Remove water soluble dirt

I Use only clean, deionised water


a Simple method: breathe upon the surface

II Use a cotton stick

 $\boldsymbol{a}$  Clean from the centre toward the edge in

concentric circles

III Inspect result



### Remove solvent soluble dirt

I Use ultra pure solvents

a Ethanol PA

**b** Acetone PA (be careful, may harm plastics and objective labels)

II Use a cotton stick, dental cotton pad

- a Cotton stick and cotton pad should only be moist, not dropping wet
- **b** Clean from the centre toward the edge in concentric circles

III Inspect result

IV Still smudgy?

a Combine the methods: breathe upon the surface and clean with cotton stick and acetone

Leica Optics Cleaning Kit
Order number **11505508**