



Living up to Life

Leica
MICROSYSTEMS

Material
Research





Life Science
Research

Application Note


Ion Beam Polishing of sample surfaces Sample Preparation for SEM

related instrument Leica EM RES102

Medical
Research



Industrial
Manufacturing



Natural
Resources



Polishing of sample surfaces

BASIC PRINCIPLE

Ion milling can be used to reduce the roughness of sample surfaces. Small angles less than 6° with respect to the sample surface are necessary. The high voltage depends on the material to be prepared.

The reason for the levelling effect is the different milling angle of flat and rough surface areas. The milling rate is lower for small angles. The rough surface area will be faster milled. The result is a flat surface.



Fig. 42: Polishing effect in case of rough surface areas

Ion polishing is often the final step of sample preparation. The prerequisite is a perfect mechanical pre-preparation as samples with deep surface scratches cannot be ion polished.

Soft materials usually have a smeared sample surface after mechanical polishing. It is necessary to remove this smeared material before ion polishing. Otherwise the above mentioned polishing effect does not work.

APPLICATION EXAMPLES

- Oil shale
- Nickel based super alloy

POLOSHING OF AN OIL SHALE SAMPLE

Purpose

Mechanical polishing of oil shale is difficult because of the different hardness of the organic and inorganic components. The surface is often not completely flat.



A final ion polishing after mechanical polishing can greatly enhance the surface quality. On the other hand, the sample surface should be as good as possible before ion polishing. Thus the combination of both steps should be the right way to achieve a good result.

PREPARATION PARAMETERS

Acceleration voltage:	6 kV
Gun current:	2.2 mA
Sample movement:	Area mode
Milling angle:	3°
Milling time:	6 h

RESULTS

The images show an oil shale sample after mechanical polishing and after additional ion polishing. Due to the different hardness of the components of the oil shale the surface quality is not good enough for imaging. An additional ion polishing of an area of about 20 mm in diameter leads to a much better surface quality. Now it is possible to observe structural details.






Fig. 1: Oil shale after mechanical polishing (upper left image) and after additional ion polishing (top right and lower image) (6 kV, 3°, 6h, and area mode)

POLISHING OF NICKEL BASED SUPERALLOY

Purpose

Usually sample surface is not completely free of damage after mechanical polishing. The depth of surface damage is three times than the grain size of the lapping foil. Ion polishing is a good tool to remove such damage. It can be used for SEM and TEM samples.

PREPARATION PARAMETERS

Acceleration voltage:	7 kV
Gun current:	2.6 mA
Sample movement:	Rotation
Milling angle:	5°
Milling time:	3 h

RESULTS

Fig. 2 shows a Nickel based superalloy after mechanical polishing and after additional ion polishing. The sample surface has a lot of fine scratches after mechanical polishing. It was possible to remove all these scratches and to achieve a nice surface after an additional 3 hours ion polishing.




Fig. 2: Nickel based super alloy (CSMX-4) after mechanical polishing (left) and after additional ion milling (7 kV, 5°, 3h, and rotation) (right)

RELATED PRODUCTS



Leica EM RES102

The statement by Ernst Leitz in 1907, “**With the User, For the User,**” describes the fruitful collaboration with end users and driving force of innovation at Leica Microsystems. We have developed five brand values to live up to this tradition: Pioneering, High-end Quality, Team Spirit, Dedication to Science, and Continuous Improvement. For us, living up to these values means: **Living up to Life.**

LIFE SCIENCE DIVISION - NANO TECHNOLOGY LNT

The Leica Microsystems Nano Technology Division's focus is to provide the most comprehensive product portfolio for the preparation of biological, medical and industrial samples for investigation in the Electron and Light Microscope. Excellent Sample Preparation is a prerequisite for perfect microscopy. **Your image starts here!**

Leica Microsystems – an international company with a strong network of worldwide customer services:

Active worldwide	Tel.	Fax
Australia · North Ryde	+61 2 8870 3500	2 9878 1055
Austria · Vienna	+43 1 486 80 50 0	1 486 80 50 30
Belgium · Diegem	+32 2 790 98 50	2 790 98 68
Canada · Concord/Ontario	+1 800 248 0123	847 405 0164
Denmark · Ballerup	+45 4454 0101	4454 0111
France · Nanterre Cedex	+33 811 000 664	1 56 05 23 23
Germany · Wetzlar	+49 64 41 29 40 00	64 41 29 41 55
Italy · Milan	+39 02 574 861	02 574 03392
Japan · Tokyo	+81 3 5421 2800	3 5421 2896
Korea · Seoul	+82 2 514 65 43	2 514 65 48
Netherlands · Rijswijk	+31 70 4132 100	70 4132 109
People's Rep. of China · Hong Kong	+852 2564 6699	2564 4163
· Shanghai	+86 21 6387 6606	21 6387 6698
Portugal · Lisbon	+351 21 388 9112	21 385 4668
Singapore	+65 6550 5999	6773 0628
Spain · Barcelona	+34 93 494 95 30	93 494 95 32
Sweden · Kista	+46 8 625 45 45	8 625 45 10
Switzerland · Heerbrugg	+41 71 726 34 34	71 726 34 44
United Kingdom · Milton Keynes	+44 800 298 2344	1908 246312
USA · Buffalo Grove/Illinois	+1 800 248 0123	847 405 0164