

White Paper

INTRAOPERATIVE
OPTICAL COHERENCE TOMOGRAPHY:
CLINICAL EXPERIENCE OF RETINAL PROCEDURES

AUTHORS:

André Maia
Ricardo Okada Nakaghi
Gabriel Costa de Andrade
Francisco Rosa Stefanini

Hermano Assis
Caio Franco da Silveira
Rafael Reis Pereira
Lucas Mares Guia Ribeiro

Table of Contents

I. Introduction.....	3
II. Case Collection.....	4
1. Macular hole repair in highly myopic patients.....	4
2. Macular hole repair and the inverted ILM flap technique.....	5
3. Treatment of retinal detachment with subretinal fibrosis.....	6
4. Tissue plasminogen activator (TPA) injection for subretinal hemorrhage.....	8
5. Autologous retinal pigment epithelium (RPE) transplantation	10
III. Surgical Findings.....	12
IV. Discussion.....	13
V. Limitations of Former Intraop. OCT Generations.....	14
VI. Conclusion.....	14
VII. References.....	15

The statements of the healthcare professionals made in this white paper reflect only their personal opinions and experiences. These statements do not necessarily reflect the opinion of any institution with whom they are affiliated.

I. Introduction

Intraoperative optical coherence tomography (intraoperative OCT) is a microscope-integrated system designed to assist ophthalmic surgeons by providing helpful information for diagnosis and surgical procedures. This valuable technology can be used in treating a variety of ophthalmic diseases and has been increasingly used in the operating room (OR).¹⁻³

Retina surgery has tremendously evolved with emergence of digital 3D heads-up viewing systems. These technologies have led to improvement in image quality, depth of focus, field of view, and also important ergonomic benefits for retina surgeons compared to the standard operating approach of viewing the procedure through the binoculars of a surgical microscope.⁴⁻⁶

The EnFocus intraoperative OCT system is built into the Proveo 8 ophthalmic microscope (Leica Microsystems). It is a high-resolution device which can be used during both anterior and posterior segment surgeries. The Leica OCT solution can be used with retina wide-angle viewing systems such as the RUV800 (Leica Microsystems) and BIOM (Möller-Wedel, Germany).

Real-time acquired OCT scans can be viewed in video mode or frame by frame via the monitor that is part of the EnFocus OCT cart solution. Also, still images from scans can be made.

In addition, the surgeon can also use this system to visualize the OCT image directly on the microscope's monitor, as well as on external monitors in the OR, thanks to the latest EnFocus OCT solution from Leica fully built into the Proveo 8.

The surgeries described in this document were performed with the 3D digitally assisted system (NGENUITY® Alcon). The images were displayed on a single screen using a HDMI standard cable connection. Thus, the surgeon was able to analyze the OCT images using a picture-in-picture view, which provided larger images as the screen could be split in two. When used together, intraoperative OCT and heads-up 3D visualization systems can improve the surgeon's decision-making regarding optimal surgical strategies aiming to improve surgical outcomes and, consequently also improving outcomes for patients.

In recent years, intraoperative OCT has become an integral part of the surgical workflow in eye surgery and herein, we share our experience of using this technology in our daily cases. To give an overview of the most interesting cases that we have performed with the EnFocus intraoperative OCT built into the Proveo 8 microscope (Leica Microsystems) we have selected cases of macular hole (MH) repair, subretinal fibrosis and submacular hemorrhage treatment and autologous retinal pigment epithelium (RPE) transplantation.




Fig. 1 Intraoperative 2D photograph of the NGENUITY® monitor showing an OCT image (EnFocus, Leica Microsystems) captured during the procedure. This image has been captured in the OR as a screen capture from one of the monitors.

II. Case Collection

1. Macular hole repair in highly myopic patients

The surgical success rate of repair of idiopathic macular holes (MH) smaller than 250 μm is about 98%, while for larger holes (400 μm or more), the closure rate is over 90%.²¹ In contrast, when highly myopic eyes are involved, MH with retinal detachment has a low surgery closure rate compared to idiopathic MH. Studies have shown considerable variations in closure rates, ranging from 60-100%, while some patients require more than one procedure to achieve closure of the MH.^{6-10,14,15}

For patients with retinal detachment and MH, the decision concerning how to treat the remaining subretinal fluid is not always clear. Specifically, there may be some uncertainty around removing the remaining subretinal fluid, and whether it should be aspirated or left. If aspiration of subretinal fluid is chosen, there may also be uncertainty around the way it should be performed – through the macular hole or a retinotomy. In addition, the aspiration procedure also has a certain risk. As the surgeon moves the cannula too close to the fluid, there is a risk of iatrogenic damage to the edges of the

macular hole, or even to the RPE. With the real-time intraoperative OCT visualization, surgeons can minimize the risk of damaging these structures during aspiration of the subretinal fluid.

The necessity and the duration of postoperative face-down positioning are still under debate.¹¹ Lorusso et al. affirmed that intraoperative confirmation of macular hole closure with intraoperative OCT could serve as a potentially useful tool for prescribing less than 24h of face-down positioning.¹²

Fig. 2 depicts the case of a 52-year-old Asian female patient, presenting with high myopia, a macular hole, and macular detachment, who underwent a 23-gauge vitrectomy and internal limiting membrane (ILM) peeling. Using intraoperative OCT, we were able to visualize the presence of subretinal fluid at the beginning of the procedure. After complete ILM peeling and fluid-air exchange, persistent subretinal fluid was identified by OCT. Using an aspiration 25-gauge cannula (Grieshaber ®), the fluid was completely drained from the subretinal space (Fig. 3). As such intraoperative OCT helped to drain the subretinal fluid, while avoiding touching the retina and RPE.




Fig. 2 Macular hole in a myopic patient with macular detachment. The intraoperative OCT image is shown on the monitor on the top right, demonstrating macular detachment.




Fig. 3: The same patient as depicted in Fig. 2. The ILM peeling and fluid-air exchange were performed. Under OCT visualization, an aspiration cannula was used to drain the submacular fluid from within the macular hole, without touching the retina and retinal pigment epithelium (RPE). The white arrow shows that the borders of the hole are approximated, and there is no macular detachment.

2. Macular hole repair and the inverted ILM flap technique

Recently, a study has demonstrated that the inverted ILM flap technique was superior in comparison with ILM peeling.¹⁶ However, visualization of the ILM flap is difficult. Intraoperative visualization of the inverted ILM flap was made possible by OCT, which we use in our daily cases. We used it in our daily cases. Fig. 4 depicts a high myopia patient being operated with inverted ILM, while intraoperative OCT shows that the macular hole is nearly closed. Also, the ILM flap is well positioned above the fovea. Fig. 5 shows the same patient depicted in Fig. 2 and 3 with another case of inverted ILM flap technique. These images show that intraoperative OCT can be used to check the ILM positioning at the end of the procedure to help surgeons achieve improved results.




Fig. 4 The intraoperative OCT is showing in detail that macular hole is nearly closed and ILM flap is in the right place.

The benefits of using intraop. OCT in macular hole repair: Key points

- > Intraoperative analysis of epiretinal membrane (to choose the optimal location to begin the peel)
- > Reduced risk of iatrogenic damage in subretinal fluid aspiration during macular detachment repair (helping to avoid touching the RPE and the retina margin)
- > Checking of the final positioning of the ILM flap (inverted ILM flap technique).
- > Evaluation of macular hole closure (serving as a potential postoperative positioning guide)




Fig. 5 The same patient as depicted in Figs. 2 and 3 with high myopia-associated macular hole and retinal detachment. The white arrow shows an inverted ILM flap placed above the macular hole.

3. Retinal detachment with subretinal fibrosis

Retinal detachment complicated with proliferative vitreoretinopathy (PVR) or subretinal fibrosis has a high rate of recurrence. The complete elimination of retinal traction is critical to ensure the success of the surgery. However, in some cases it is challenging to establish whether the PVR has been removed. Moreover, it might be difficult to establish differences between PVR, and the subretinal fibrosis band. Using intraoperative OCT can offer important advantages to surgeons during these procedures. Fig. 6 shows a case of chronic inferior retinal detachment.




Fig. 6 After removing the PVR, the inferior retina was detached, revealing the subretinal fibrosis (white arrow). The EnFocus intraoperative OCT image shows an area with some subretinal fluid in the inferior edge of the retina (depicted by the arrowhead on the top left).

The PVR was removed and the perfluorocarbon liquid (PFCL) was injected into the vitreous cavity. It was possible to observe that the inferior retina had a shallow detachment. Moreover, the intraoperative OCT had confirmed that the inferior retina was not flat and showed a subretinal fibrotic band.

The PFCL was removed and an endodiathermy probe and 25-gauge forceps were used to remove the fibrosis (Fig. 7), flattening the retina. A fluid-air exchange was then performed, and a gas tamponade was used to attach the retina (Fig. 8). The use of intraoperative OCT contributed to the visualization of the subretinal band, its removal and further observation of retinal reattachment. In this case, intraoperative OCT was vital for distinguishing subretinal fibrosis from PVR and ensuring retinal reattachment.




Fig. 7 An endodiathermy probe was used for retinotomy and subretinal fibrosis was removed with 23-gauge forceps placed above the macular hole.




Fig. 8 The same patient as depicted in Figs. 5 and 6. Attached retina after removing subretinal fibrosis band. EnFocus intraoperative OCT confirmed the retinotomy carried out using an endodiathermy probe (depicted by the arrow).




Fig. 9: A chronic retinal detachment with a subretinal fibrotic band (depicted by the white arrow). After core vitrectomy and posterior hyaloid detachment, the retina still hasn't completely settled and subretinal bands are observed. EnFocus intraoperative OCT image at top right of the monitor shows the subretinal fibrotic band (depicted by the arrowhead). Subretinal fibrosis was removed with 23-gauge forceps.

4. Subretinal hemorrhage

Subretinal surgeries have been increasingly performed in recent years, and the injection of tissue plasminogen activator (TPA) can be facilitated by intraoperative OCT. Since the treatment for genetic conditions continues to evolve, mastering the technique of subretinal injections is crucial for the retinal surgeon.

Extensive subretinal hemorrhages may be caused by any macular condition associated with choroidal neovascularization. The most common cause is age-related macular degeneration (AMD). Other cases have been described as polypoidal choroidal vasculopathy (PCV), pathologic myopia, angioid streaks, trauma, and retinal arterial macroaneurysm.

Due to the prognosis of submacular hemorrhages, many different surgery modalities have been used in an attempt to improve visual

outcomes. The result of a retrospective study suggests that pars plana vitrectomy, subretinal TPA with pneumatic displacement lead to better visual outcomes compared to the natural course of this disease.¹⁷

In our clinic, a case of extensive submacular hemorrhage was managed using assisted OCT (Fig. 10). Specifically, a 23-gauge pars plana vitrectomy with posterior vitreous detachment and ILM peeling was performed. A 36-gauge cannula connected to a fluid injector was inserted near the temporal vascular arcades and BSS fluid was injected to create a subretinal space. The retinal posterior pole was detached. After that, subretinal TPA with anti-VEGF was injected with a 36-gauge cannula (25ug/0,1 ml). Fluid-air exchange was performed, pneumatic displacement was then realized with C3F8 15% and the patient was instructed to maintain the head in an upright position.

In addition, intraoperative OCT was used to observe the interaction between the TPA cannula and the retina. The subretinal positioning of the cannula was confirmed with intraoperative OCT; the fluid was safely injected, and the macula was detached (Fig. 10). At this point, a possible complication can arise during insertion of the cannula into the sub-choroidal space, as we can see in the example shown below in Fig. 11.

The wrong maneuver was promptly identified and intraoperative OCT was used to correct it (Fig. 11). Another potential complication arising after submacular injection is creation of an iatrogenic macula hole. However, OCT use allows preventing misdiagnosis of this complication.



Fig. 10 Subretinal injection with a 36-gauge cannula. First, a balanced salt solution (BSS) was injected to detach the macular region, then, an anti-VEGF solution with TPA was injected in the same location. Finally, fluid air exchange with a 23-gauge cannula was performed and pneumatic displacement was carried out.




Fig. 11 Iatrogenic choroidal detachment visualized using OCT (depicted by white arrows). The cannula was then removed and a choroidal misdirection of TPA was avoided thanks to the EnFocus intraoperative OCT visualization.

5. Autologous RPE transplantation

RPE transplantation has been investigated as an option for the treatment of retinal damage in exudative or atrophic maculopathies such as RPE tears in age-related macular degeneration (AMD), massive subretinal hemorrhage, and geographic atrophy.¹⁸ It is an advanced technical procedure that requires accurate surgical dexterity. The additional use of intraoperative OCT can offer some benefits for RPE transplantation, since this technology provides real-time feedback about the correct positioning of the patch.

Figs. 11 to 14 depict a RPE transplantation procedure performed in an 86-year-old woman with choroidal neovascularization related to AMD. During the procedure, intraoperative OCT helped to confirm the precise positioning of the RPE transplant. Also, intraoperative OCT showed a PFCL bubble in the inferior border of the patch, with liquid removal also guided by this technology. Thus, intraoperative OCT provided real-time feedback for the surgeon, preventing further complications.




Fig. 12 After complete removal of vitreous core and base with 23-gauge vitrectomy, ILM peeling with 23-gauge forceps was performed. The temporal retina was folded onto the nasal retina and the RPE patch was created using endodiathermy. The EnFocus intraoperative OCT image scan shows the RPE and choriocapillaris on the autologous patch.




Fig. 13 Perfluorocarbon liquid (PFCL) was used for creating subretinal space and an RPE patch was excised and moved toward the macular region. The figure shows the RPE patch being placed in the macular region. EnFocus intraoperative OCT was used to scan the posterior pole to obtain the position of the transplant.




Fig. 14 The intraoperative OCT scan performed with EnFocus revealed a PFCL bubble on the inferior edge at the margin of the patch. Since the presence of the PFCL bubble and hemorrhage are risk factors for failure of the procedure, we decided to remove the bubble. The PFCL above the retina was aspirated: the temporal retina was folded again and the subretinal PFCL at the border of the patch was then aspirated. Subsequently, the retina was replaced using PFCL.




Fig. 15 The well positioned autologous RPE transplant is shown on the Ngenuity® monitor.

III. Surgical Findings

EnFocus intraoperative OCT can help to reveal important details to help surgeons during procedures. Some of these interesting features shown by intraoperative OCT have caught our attention. As we expected, forceps and probes are opaque and do not allow image formation. However, the light beam can pass through the brilliant blue dye (Membraneblue Dual® DORC) without any restrictions. As a result, the scattering of the beam can produce clear retina images (Fig. 16).

We also observed another interesting feature during the fluid-air exchange. The intraoperative OCT image could be produced within the air, helping us determine the presence of subretinal fluid underneath (Fig. 17).




Fig. 16 Brilliant blue dye is used worldwide for staining internal limiting membranes. The EnFocus intraoperative OCT image is depicted in the picture-in-picture. The OCT light beam can pass through the dye with no interference.




Fig. 17 Rhegmatogenous retinal detachment with retinal macroaneurysm. After the fluid-air exchange, it was possible to observe remaining subretinal fluid on the macular area using EnFocus intraoperative OCT. Moreover, the OCT image was produced despite the cavity being filled with air.

IV. Discussion

The usefulness of intraoperative OCT in ophthalmic procedures has been evaluated in numerous studies. Its positive impact on surgical decision-making has been reported in DISCOVER and PIONEER studies (involving prospective, consecutive case series).

The PIONEER study demonstrated the feasibility of use of intraoperative OCT imaging for numerous anterior and posterior segment ophthalmic procedures in the operating room²⁰

In this study, a portable spectral domain OCT (Biotigen SDOIS, Research Triangle Park, North Carolina, USA) was stabilized with a microscope-mount system. Since the OCT was not integrated into the microscope and the surgery needed to be paused to acquire images, the technology used in the PIONEER study was unable to visualize the interactions between the tissue and instruments used.

The DISCOVER study evaluated the utility of intraoperative OCT and assessed surgical outcomes in patients with retinal detachment (RD).¹⁹ The study involved the use of the RESCAN 700 microscope-integrated OCT (Carl Zeiss Meditec, Inc, Oberkochen, Germany), EnFocus intraoperative OCT (Biotigen/Leica Microsystems, Wetzlar, Germany) and Cole Eye Institute's microscope-integrated intraoperative OCT (Cleveland Clinic Foundation, Cleveland, Ohio). In terms of RD repairs, the study points out the potential benefits of OCT in treatment of complex pathologies, such as PVR, panuveitis, giant retinal tears and recurrent retinal detachment. Specifically, intraoperative OCT was demonstrated to provide useful feedback in 50% of complex cases, and in 22% of regular cases (non-complicated RD). The study also showed that in 12% of cases, intraoperative OCT data directly altered surgical decision making. Despite the lack of a control group, the study showed a reoperation rate of 25% in complex cases, which closely corresponds the overall reported ranges of 13-30%.¹⁹

In our service, intraoperative OCT was used in a variety of surgical cases and allowed us to discover new insights in retinal surgery. In addition, the use of intraoperative OCT also provided a better understanding of retina structures over the course of the surgery.

For example, retinal areas under suspicion of containing subretinal fluid after fluid-air exchange could be reviewed, helping to improve the surgeon's confidence during the procedure. In addition, intraoperative OCT was also able to help the surgeon ensure the correct positioning of the inverted ILM flap, since it is not possible to detect the integrity of the ILM flap intraoperatively through any other means. Another example of the benefits of intraoperative OCT use is improved intraoperative visualization of macular hole closure. A recently published study has demonstrated that intraoperative OCT-based identification of macular hole closure could serve as a safe and helpful tool when prescribing short-term postoperative prone posturing.¹²

The combination of intraoperative OCT with a 3D-assisted display system has also brought forth new possibilities for ophthalmic surgery. Through a simple connection with an HDMI cable, the image generated by EnFocus can be projected (using picture-in-picture) on the same 3D screen used by the surgeon during surgery. This makes it possible to maintain optimal ergonomics (heads-up), while reducing intraoperative time to analyze the images (since there is no need to look at another screen). Moreover, this technology allows to increase the size of the image (splitting the screen/zoom in and out) and assisting the surgeon in the decision-making process.

V. Limitations of Former Intraop. OCT Generations

The integration of previous generations of intraoperative OCT technology into the surgical workflow involved some challenges, and the help of an imaging assistant was often needed to assist the surgeon in handling the intraoperative OCT modality during surgery. For example, surgeons required assistance during procedures when it was necessary to change a scan pattern or optimize image quality of live OCT images.

Despite the possibility of using a preprogrammed footswitch to independently control movements in the X-Y-Z axis, capture still images, and acquire scans, surgeons relied on their assistants during the surgery to support them in efficient operation of the OCT.

Recent advances in the integration of OCT into the surgical workflow have succeeded in closing this gap. The next-generation EnFocus intraoperative OCT built-into the Proveo 8 surgical microscope (Leica Microsystems) provides the surgeon with

VI. Conclusion

Retinal surgery continues to evolve, while surgeons continuously strive to further optimize and develop their surgical skills and techniques to improve clinical outcomes. With improved integration of intraoperative OCT modalities into the surgical workflow, the microscopic view becomes supplemented with bright and sharp images of previously hidden subsurface details, which enable surgeons to better understand ocular pathology. This novel technology will further transform retina surgery. Furthermore, intraoperative imaging in real-time provides a tremendous advantage for surgeons, as it allows them to observe the way retinal tissue reacts to surgical maneuvers. Moreover, intraoperative OCT provides surgeons with high-resolution images and immediate confirmation of surgical maneuvers, allowing observation of tissue changes at every step of the procedure. This allows surgeons to gain confidence and improve patient outcomes.

The combination of intraoperative OCT imaging with a 3D visualization system brings further ergonomic benefits to the

maximum freedom. Using this cutting-edge technology, surgeons can work smoothly and independently. In addition, there is no more need to invest time and effort in optimizing and aligning OCT imaging during live procedures. The new automatic image optimization and tracking of the EnFocus OCT help to assure that surgeons will always have consistent, optimized OCT imaging available whenever it is needed.

Furthermore, using a technician to operate OCT during procedures is not required anymore, as surgeons can switch from the microscope view to an OCT view effortlessly, with a tap on the footswitch or on the handgrip. The full integration of OCT into the Proveo 8 microscope also facilitates recording of procedures, which can be controlled with the microscope footswitch, handgrip or even via touchscreen control of the microscope's monitor. The latter can be used to facilitate touch-friendly gestures on the user-interface screen, such as manual scrolling and using the two-finger pinch to resize images, as well as pointing to specific locations, moving and rotating OCT scans.

surgical workflow, while minimizing interruptions and capturing quality OCT images. In our surgical practice, we display both

images together on a single monitor. Moreover, intraoperative OCT is especially useful in cases with pre-retinal inflammatory precipitates (PIP), as this technology visually aids in performing sensitive maneuvers that require a high level of precision, while helping to avoid damaging retinal tissue.

The use of intraoperative OCT during complex eye surgery, and especially in retina surgery, opens up new possibilities for improved clinical treatment and potentially improving postoperative results.

The clinical utility of intraoperative OCT technology will become more evident as more prospective clinical trials are conducted, which will provide additional evidence of its efficacy and the associated benefits for patient outcomes.

VII. References

1. Ehlers JP, Goshe J, Dupps WJ, et al. Determination of feasibility and utility of microscope-integrated optical coherence tomography during ophthalmic surgery: The DISCOVER study RESCAN results. *JAMA Ophthalmol.* 2015;133(10):1124-1132
2. Ehlers JP, Modi YS, Pecun PE, et al. The DISCOVER Study 3-Year Results: Feasibility and Usefulness of Microscope-Integrated Intraoperative OCT during Ophthalmic Surgery. *Ophthalmol.* 2018;125(7):1014-1027
3. Hattenbach LO, Framme C, Junker B, Pielen A, Agostini H, Maier M. Intraoperative real-time OCT in macular surgery. *Ophthalmologie.* 2016 Aug;113(8):656-62
4. Eckardt C, Paulo EB. Heads-up surgery for vitreoretinal procedures: An experimental and clinical study. *Retina.* 2016. 36:137-147
5. Weinstock RJ, Ainslie-Garcia MH, Ferko NC, Qadeer RA, Morris LP, Cheng H, Ehlers JP. Comparative Assessment of Ergonomic Experience with Heads-Up Display and Conventional Surgical Microscope in the Operating Room. *Clin Ophthalmol.* 2021;15:347-356.
6. Palácios RM, de Carvalho ACM, Maia M, Caiado RR, Camilo DAG, Farah ME. An experimental and clinical study on the initial experiences of Brazilian vitreoretinal surgeons with heads-up surgery. *Graefes Arch Clin Exp Ophthalmol.* 2019 Mar;257(3):473-483
7. Kuriyama S, Matsumura M, Harada T, Ishigooka H, Ogino N. Surgical techniques and reattachment rates in retinal detachment due to macular hole. *Arch Ophthalmol* 1990; 108:1559–1561.
8. Kelly NE, Wendel RT. Vitreous surgery for idiopathic macular holes. Results of a pilot study. *Arch Ophthalmol* 1991;109:654–659.
9. Akiba J, Konno S, Yoshida A. Retinal detachment associated with a macular hole in severely myopic eyes. *Am J Ophthalmol.* 1999;128:654–655
10. Ikuno Y, Sayanagi K, Oshima T, et al. Optical coherence tomographic findings of macular holes and retinal detachment after vitrectomy in highly myopic eyes. *Am J Ophthalmol* 2003;136:477–481.
11. Chow DR, Chaudhary KM. Optical coherence tomography-based positioning regimen for macular hole surgery. *Retina.* 2015. 35:899-907.
12. Lorusso M, Ferrari LM, Cicinelli MV, et al. Feasibility and Safety of Intraoperative Optical Coherence Tomography-Guided Short-Term Posturing Prescription after Macular Hole Surgery. *Ophthalmic Res.* 2019. DOI: 10.1159/000501561
13. De Giacinto C, Pastore MR, Cirigliano G, Tognetto D. Macular Hole in Myopic Eyes: A Narrative Review of the Current Surgical Techniques. *J Ophthalmol.* 2019 Mar 11;2019:3230695
14. Frisina R, Gius I, Palmieri M, Finzi A, Tozzi L, Parolini B. Myopic Traction Maculopathy: Diagnostic and Management Strategies. *Clin Ophthalmol.* 2020 Nov 2;14:3699-3708
15. Matsumae H, Morizane Y, Yamane S, et al. Inverted Internal Limiting Membrane Flap versus Internal Limiting Membrane Peeling for Macular Hole Retinal Detachment in High Myopia. 2020. *Ophthalmol.*, Vol 127, Issue 10, October 2020, Pages 1282-1286.
16. Mete M., Alfano A., Guerriero M., et al. Inverted internal limiting membrane flap technique versus complete internal limiting membrane removal in myopic macular hole surgery. *Retina.* 2017;37(10):1923–1930
17. Juncal VR, Hanout M, Altomare F, et al. Surgical management of submacular hemorrhage: experience at an academic Canadian centre. 2018. *Can J Ophthalmol.* Aug;53(4):408-414
18. Parolini B, Salvatore A, Pinackatt SJ, et al. Long-results of autologous retinal pigment epithelium and choroid transplantation for the treatment of exudative and atrophic maculopathies. 2020. *Retina.* 40:507-520.
19. Abraham JR, Srivastava SK, Le TK, et al. Intraoperative OCT-Assisted Retinal Detachment Repair in the DISCOVER Study: Impact and Outcomes. *Ophthalmol Retina.* 2020 Apr;4(4):378-383
20. Ehlers JP, Dupps WJ, Kaiser PK, et al. The Prospective Intraoperative and Perioperative Ophthalmic ImagiNg With Optical CoherEncE TomogRaphy (PIONEER) Study: 2-Year Results. *Am. J. Ophthalmol.* 2014.Vol 158 (5); Nov; 999-1007.
21. Gaudric A, Tadayoni R. (2019) "Macular hole", in Schachat AP (ed) *Ryan's Retina.* Elsevier, p. 2219



FOCUS ON PERFECTION

EnFocus intraoperative OCT built into the Proveo 8 microscope provides

- > Greater Insight
- > Immediate Confirmation
- > Maximum Freedom

MC-0003034 - 04.08.2021 - EN - Copyright © 2021 Leica Microsystems (Switzerland) AG. All rights reserved. Subject to modifications. LEICA and the Leica Logo are registered trademarks of Leica Microsystems IR GmbH. Other trademarks shown herein are the property of their respective owners.



Proveo 8 is a Class I surgical microscope




Leica Microsystems (Schweiz) AG
Max Schmidheiny-Strasse 201
9435 Heerbrugg, Switzerland



EnFocus OCT is a class IIa medical device



Leica Microsystems NC, Inc.
4222 Emperor Blvd, Suite 390,
Durham, NC 27703, USA



Leica Microsystems CMS GmbH
Ernst-Leitz-Strasse 17-37
35578 Wetzlar, Germany

Not all products or services are approved or offered in every market. Approved labeling and instructions may vary between countries. Please contact your local Leica representative for details.

CONNECT WITH US!



Leica Microsystems (Schweiz) AG · Max Schmidheiny Strasse 201 · CH-9435 Heerbrugg

T +41 71 726 3333 · F +41 71 726 3399

www.leica-microsystems.com