

Leica HCS LSI

High Content Screening Automation

Technical Documentation



Leica HCS LSI - High Content Screening Automation

Explore the new dimensions of imaging!

High content is provided by true confocal high resolution point scanner technology for pre-screen and secondary screens. Fast analysis is performed by digital camera imaging. Maximum flexibility is generated by the adaptive zoom technology, combining the innovative optical zoom with confocal zoom for universal applications. Large working distance macro objectives as well as micro objectives offer perfect system adaptation to your experiments.

The High Content Screening Automation software enables efficient screening and easy automation. Computer Aided Microscopy allows external system control and turns your Leica imaging system into an intelligent microscope.


For image analysis, the OME.TIF format provides highest compatibility with existing solutions. The Leica Edition of Definiens Developer XD even enables 4D image analysis.



Leica HCS LSIApplication Solution



Automation: Leica LAS AF Matrix M3



Analysis: Leica Edition Definiens Developer XD

Specifications

True spectral confocal poi	int scanning system		
Scanner	Method		True confocal point scanner
	Confocal channels		1
	Scanner		Galvo, [x, y]
	Sequential scan		yes
	Channels		1- 8, sequential multiplexing
.aser	Laser type		solid state
	Laser		max 4
	Laser excitation wavelength	[nm]	405, 488, 561, 635
	Excitation attenuation		AOTF, direct modulation
	Excitation attenuation control		automated
	Range	[%]	0-100
pectral detection	Spectral detection		yes
	Band selection		continuously variable
	Bandwidth	[nm]	430 – 750
	Spectral resolution	[nm]	5 nm
	Detector		1, direct coupling
	Detector type		ultra high dynamic PMT
ransmitted light	Transmitted light detector		yes
	Detector type		ultra high dynamic PMT
Resolution	Range (min – max)	[pixel²]	128 – 2048
	Scan formats	[pixel]	128, 256, 512, 1024, 2048
	Image bit depth	[bit]	8 or 12, switchable
Beam splitter	Туре	, ,	high performance dichroics
	Beam splitter 1, wavelength	[nm]	405, 488, 561, 635
	Beam splitter 2, wavelength	[nm]	405, 588, 561
	Beam splitter 3, ND	[%]	30/70
	Beam splitter exchange		automated
inhole	Pinhole type		motorized, variable
	Range (min - max)	[µm]	35 - 600
	Pinhole adjustment	[%]	0 – 100
	Control		automated via GUI
can modes			
	2D	Line, time, area	xt, xy, xz
	3D	Volume	xyz; xzy
		Area, time	xyt; xzt
		Area, lambda	χγλ; χzλ
	4D	Volume, time	xyzt; xzyt
		Area, lambda, time	xyλt; xzλt
		Volume, lambda	хуλг
		Area, lambda, time	xyλt
	5D	Volume, lambda, time	xyzλt
peed	Speed mode	fu-1	uni-, bi-directional
	Line speed range	[Hz]	400, 600, 800, bi: 1200, 1600
	max @ 128² bi-directional	[f/s]	7.0
	standard @ 512 ²	[f/s]	2.0
	min @ 2048 ²	[f/s]	0.36
0V	Field of view (diagonal)	[mm]	16
ower supply	Power supply integration		yes
.,,	Туре		auto select
	Voltage range	[V]	100-240 AC
	Digital comore		DFC365 FX
Digital camera imaging	Digital camera Number of dichroics		4 + empty, manual turret

vvorkstation				
External computer	Processor		Intel® Xeon® Quad (Core
	Memory	[GB]	12	
	HD-Size	[GB]	2000	
	Operating system		Windows 7®	
Interfaces	USB		8	
	FireWire		4 (1 x 1394a, 3 x 139	4b)
	Ethernet		1	
	DVI, HDMI		1/1	
Monitor	Graphics resolution	[Pixel]	2 x (2560 x 1600)	
	Monitor size		2 x 19"/30" (48 cm/7	'6 cm)
Micro & macroscope stand	i			
Characteristics	Туре	upright	LSI6000, base	
	Application		Micro and Macro	
	Focusing drive		motorized	
	Head Travel Range	[nm]	150	
Illumination	Transmitted light, intensity control		automated and mar	nual
	Fluorescence illumination		EL6000	
	Contrasting		Rottermann tilted ill	umination
Workspace	Height, Depth, Width	[mm]	180, 420, 555	
	Wing door access, open	[°]	180	
Laser safety	Laser safety acrylic glass box		yes	
	Laser safe tube		yes	
Micro manipulation	Manipulator type		Eppendorf, motorize	ed
	Mounting		inside, variable	
Environmental control	Sample light protection		yes	
	Temperature		yes	
	CO,		yes	
	Humidity		yes	
	Gas cover for galvo stage		yes	
Z-drives	Туре	Drive	Travel range [µm]	Z-resolution [µm]
	SuperZ stage	Galvanometer stage	500	0.01
	Fine focus	Stepper motor	10,000	5.0
		•		

Leica HCS LSI – Optics

Motor focus

	1			
Adaptive Zoom Technology				
	Zoom types		2	
Confocal zoom	Туре		Confocal scanner in	itegration
	Magnification range		1x - 58x	
	Zoom increment		0.1	
	Zoom control		motorized, continuo	
	NA		0.3 – 1.30, objective	independent
0-4:1	T		70 400 4	71C A DO A
Optical zoom	Type	or de ble	Z6 APO A	Z16 APO A
	Magnification range	variable	0.57 – 3.6x	0.57 – 9.6x
	Zoom increment		0.01	0.01
	NA, w objective 5x	variable	0.10 - 0.50	0.09 - 0.50
	Zoom control		motorized, continuo	usly variable
	Focus control		motorized fine focus	s optics
	Optical z-positioning		motorized	
	Diaphragm		motorized	
Objectives	Magnification		NA, objective	Working distance [mm]
Macro objectives	1x		0.117	97.0
	2x		0.234	39.0
	5x		0.5	19.0
Minus abinations	10		0.0	0.00
Micro objectives⁵	10x		0.3	0.30
	20x		0.6	0.16
	40x		0.8	0.16
4	63x		1.3	0.15
4				

Stepper motor

150,000

1.0

Paran	neter		Fie	ld of v	iew [m	m]			0r	otical n	nagnifi	cation	[x]	
Zoom	type	Z1	6 APO	Α	Z	6 APO	A		Z1	16 APO	Α	Z	6 APO	A
Objec	tive	1x	2x	5x	1x	2x	5x		1x	2x	5x	1x	2x	5x
	0.6	22.0	11.0	4.4	22.0	11.0	4.4	0.6	0.7	1.4	3.5	0.7	1.4	3.5
	0.8	15.7	7.8	3.1	15.7	7.8	3.1	0.8	1.0	2.0	5.0	1.0	2.0	5.0
	1.0	12.5	6.3	2.5	12.5	6.3	2.5	1.0	1.2	2.5	6.2	1.2	2.5	6.2
8	1.3	10.0	5.0	2.0	10.0	5.0	2.0	1.3	1.6	3.1	7.8	1.6	3.1	7.8
Ling.	1.6	7.8	3.9	1.6	7.8	3.9	1.6	1.6	2.0	4.0	9.9	2.0	4.0	9.9
settings	2.0	6.3	3.1	1.3	6.3	3.1	1.3	2.0	2.5	5.0	12.4	2.5	5.0	12.4
E	2.5	5.0	2.5	1.0	5.0	2.5	1.0	2.5	3.1	6.2	15.5	3.1	6.2	15.5
zoom	3.2	3.9	2.0	8.0	3.9	2.0	8.0	3.2	4.0	7.9	19.8	4.0	7.9	19.8
- E	3.6	3.5	1.7	0.7	3.5	1.7	0.7	3.6	4.5	8.9	22.3	4.5	8.9	22.3
Optical	4.6	2.7	1.4	0.5				4.6	5.7	11.4	28.5			
0	5.0	2.5	1.3	0.5				5.0	6.2	12.4	31.0			
	6.3	2.0	1.0	0.4				6.3	7.8	15.6	39.1			
	8.0	1.6	0.8	0.3				8.0	9.9	19.8	49.6			
	9.2	1.4	0.7	0.3				9.2	11.4	22.9	57.2			

Note: The dark green fields mark the recommended range for 3D imaging at confocal zoom 1. For Z16 APO A, NA increases up to 6.3x optical zoom.

Magnification range, optical zoom

max. 0.6x - 9.2x

Macro objectives with	Combination of optical and confocal zoom	yes
optical and confocal zoom	Magnification range of optical and confocal zoom, total	max. 0.6x - 533.6x

iviicro objectives	iviagnification range of confocal zoom, total						
with confocal zoom	Parameter Field of view [mm]						
	Adapter 1x	15 6904 623					
	Objective 10x 20x 40x 63x						
	Confocal zoom 1	1 60	0.80	0.40	0.25		

Optical parameters
Macro objectives
with optical zoom

IX — 58X							
Optical magnification [x]							
	15 6904 623						
	10x	20x	40x	63x			
	10x	20x	40x	63x			

Leica HCS A – High Content Screening Automation

Leica HCS A platform info			
	System support ¹		TCS LSI, HCS LSI
	LAS AF version		2.5.0 or higher
Image acquisition	Imaging technologies		True confocal point scanner
-			Digital imaging
	Supported cameras		DFC365 FX
Multicolor	Confocal	No. of colors	8
	Camera	No. of colors	1
Transmitted light	Confocal		yes, optional
·	Camera		yes, optional
Motorized stage	Scanning stage	15 6905 202	Included in HCS LSI only
·	xy-travel range ²		127 x 83 mm
Export formats	Image types ³		TIF, OME.TIF, LIF
_	Image data format		OS platform independent
Network	Protocol		TCP/IP
	Administration		Local system admin
Remote system control	Control via network		Yes, with CAM
,	Control interface		Computer Aided Microscopy, CAM
System requirements	Platform		TCS LSI, HCS LSI
	Operating system		Windows XP®, SP3, Windows 7®
	Recommended		
	Processor speed	[MHz]	3
	Memory	[GB]	4
	Hard disk	[GB]	500
	Network		yes
Limitations		Hardware excluded	TCS SP5 X, TCS SP5 MP, TCS STED, DMI6000 CFS
			TCS SMD FCS, FLIM, FLCS, WLL
		Software excluded	FRAP, FRET, Electrophysiology
			. , , , , , , , , , , , , , , , , , , ,

Leica HCS A software			
Imaging automation	Licenses		Included in the HCS LSI product
	LAS AF MATRIX Mosaic Advanced	156602501	yes
	LAS AF MATRIX Mosaic+Multiwell Advanced	156602502	yes
	LAS AF MATRIX Mosaic Full Version	156602504	yes
	LAS AF MATRIX Multiwell Full Version	156602505	yes
	LAS AF MATRIX Full Version w/o CAM	156602511	yes
Remote system control	LAS AF MATRIX Developer Entry	156602512	upgrade option, not included
	LAS AF MATRIX Developer Full w. CAM	156602514	upgrade option, not included
Accessory tools	LAS AF MATRIX Single Object Tracking	156602507	upgrade option, not included
	LAS AF MATRIX Z-Drift Compensator	156602509	upgrade option, not included

Leica HCS LSI – System Performance

Leica HCS LSI system performance							
	Acquistion speed measureme	ent	Avg. time per well in	Avg. time per well incl. xy-travel, 96-well plate			
True confocal screen	HCS LSI w. confocal system						
Parameter	No. of channels incl. TLD	No. of colors	No. of z-sections	Format	Confocal [min/well]		
Primary screen	2	1	1	256 x 256	0.04		
	2	1	1	512 x 512	0.04		
	2	1	1	1024 x 1024	0.06		
Secondary screen	2	1	1	2048 x 2048	0.06		
	2	1	10	512 x 512	0.36		
	2	1	30	512 x 512	0.37		
	2	1	100	512 x 512	0.89		

Digital camera screen	HCS LSI w. DFC360 FX, binning: none, image depth: 12 bit					
	No. of channels, no TLD	No. of colors	No. of z-sections	Format	Camera [min/well]	
	1	1	1	1392 x 1040	0.06	
	1	1	10	1392 x 1040	0.09	
	1	1	30	1392 x 1040	0.13	
	1	1	100	1392 x 1040	0.31	

Leica HCS LSI - Compatible Image Analysis and Control


Image analysis & remote control software options	
Compatible software (examples)	Not included in the packages
Programming languages	C++, C#, VB, Lab VIEW™, MATLAB™
Image manipulation	Adobe Photoshop®
lmage analysis software⁴	ImageJ with LOCI plug-ins for OME import
	MetaMorph®, MM AF® TIFF import
Leica Edition of Definiens	Not included in the packages
Image analysis software	Leica Edition of Definiens Developer XD
Import	OME.TIF, LIF, Metadata
Programming & Plug-ins	yes
Report out	yes

Annotations

- 1 Supports technology within the range of the product specification.
- $2\,\mbox{For all}$ sample carrieres, a test is recommended.
- 3 Open Microscopy Environment (OME) is a multi-site collaborative effort among academic laboratories and a number of commercial entities that produces open tools to support data management for biological light microscopy. Designed to interact with existing commercial software, OME source code is available under GNU public copyleft licenses. OME is developed as a joint project between research-active laboratories at the Dundee, NIA Baltimore and Harvard Medical School and LOCI
- 4 C++ is a programming language standardized by ISO. C# is a programming language developed by Microsoft, Inc. Lab VIEW™ is a registered trademark of NI National Instruments Inc. MATLAB™ is a registered trademark of The MathWorks™, Inc. Adobe Photoshop® is a registered trademark of Adobe Systems® Incorporated. ImageJ is a public domain Java image processing program inspired by National Institutes of Health, NIH. Windows® XP is a registered trademark of the Microsoft® Corporation. Definiens® is a registered trademark of Definiens AG. Meta-Morph® is a registered trademark of MDS Analytical Technologies.
- 5 The software supports switching between two macro objectives. For micro objectives, an adapter is required. For micro objectives, only the confocal zoom is applicable.

Room Requirements


n .	D 1.1.				
Power supply	Power supply integration		yes		
	Type		auto select		
	Voltage range	[V]	100 – 240 AC		
	Power consumption	[VA]	800		
	Independent circuits	[no.]	1		
	Frequency	[Hz]	50/60		
	Fuse: standard	[A]	10		
	Note: The optimal optical performance can only be achieved on stable room floors. Concrete floors are required. Others, i.g. wooden floors, are not suitable.				
Environment	Humidity	[%]	10 – 80		
	Operating temperature	[°C]	18 – 30		
	Guaranteed stability	[°C]	23 +/- 2		
Load capacity and weight	Confocal unit, max.	[kg]	75		
	Microscope, max.	[kg]	45		
	System	[kg]	90		
	Static floor load	[kg/m²]	200		




Visible radiation:





LASER RADIATION VISIBLE AND INVISIBLE - CLASS 3B AVOID DIRECT EXPOSURE TO BEAM < 500mW 350-700nm IEC 60825-1: 2007



Additional Leica HCS LSI options















15 6901 303 Leica HCS LSI m. Basis CSQ RYB **15 6901 321** Leica HCS LSI m. Basis CSQ RYBV











Leica HCS LSI Product Overview

D 1 10 :			
Product Overview			
Platforms	System support		HCS, LSI, TCS LSI
Image acquisition	Technology		True spectral confocal
			Digital camera option
	Available laser lines, solid state		405, 488, 532, 561, 635
	Pre-scan and secondary scan		yes
Motorized stage	Scanning stage	15 6905 202	yes
	Multiwell xy-travel range	127 x 83 mm	yes
Environmental control	Climate chamber		
Imaging automation	LAS AF MATRIX Mosaic Advance	15 6602 501	yes
	LAS AF MATRIX Mosaic + Multiwell Advanced	15 6602 502	yes
	LAS AF MATRIX Mosaic Full Version	15 6602 504	yes
	LAS AF MATRIX Multiwell Full Version	15 6602 505	yes
	LAS AF MATRIX Full Version w/o CAM	15 6602 511	yes
Workstation	Power PC with Intel Core Duo-Processor		yes
	Operating system	Windows 7®	yes
	Monitor TFT 19" (48 cm)	2	yes
	Keyboard, Mouse	1	yes



- Confocal scanhead
- Optical zoom, motorized
- Motor focus drive
- Laser safety chamber
- 4 5 6 7 8 Wing doors
- SuperZ Galvo stage
- xy-stage control device
- Macroscope touch control
- 9 Confocal control panel
- 10 Leica EL6000 fluorescence illumination
- **①** Keyboard, mouse
- Monitor
- 13 Digital camera
- Confocal laser supply unit
- 15 Workstation
- Antivibration table, passive
- Computer table
- Workspace
- Motorized xy-stage
- 16 17 18 19 20 21 Heat pipe adapter
- Transmitted light detector
- Micro- and macroscope stand





Please note that as result of the continuous development of our products the data and final appearance can vary from the information provided in this document.

Subject to modifications. LEICA and the Leica Logo are registered trademarks of Leica Microsystems IR GmbH.

"With the user, for the user" Leica Microsystems

Leica Microsystems operates globally in four divisions, where we rank with the market leaders.

• Life Science Division

The Leica Microsystems Life Science Division supports the imaging needs of the scientific community with advanced innovation and technical expertise for the visualization, measurement, and analysis of microstructures. Our strong focus on understanding scientific applications puts Leica Microsystems' customers at the leading edge of science.

Industry Division

The Leica Microsystems Industry Division's focus is to support customers' pursuit of the highest quality end result. Leica Microsystems provide the best and most innovative imaging systems to see, measure, and analyze the microstructures in routine and research industrial applications, materials science, quality control, forensic science investigation, and educational applications.

Biosystems Division

The Leica Microsystems Biosystems Division brings histopathology labs and researchers the highest-quality, most comprehensive product range. From patient to pathologist, the range includes the ideal product for each histology step and high-productivity workflow solutions for the entire lab. With complete histology systems featuring innovative automation and Novocastra™ reagents, Leica Microsystems creates better patient care through rapid turnaround, diagnostic confidence, and close customer collaboration.

Medical Division

The Leica Microsystems Medical Division's focus is to partner with and support surgeons and their care of patients with the highest-quality, most innovative surgical microscope technology today and into the future.

The statement by Ernst Leitz in 1907, "with the user, for the user," describes the fruitful collaboration with end users and driving force of innovation at Leica Microsystems. We have developed five brand values to live up to this tradition: Pioneering, High-end Quality, Team Spirit, Dedication to Science, and Continuous Improvement. For us, living up to these values means: Living up to Life.

Active worldwide

Australia:	North Ryde	Tel. +61 2 8870 3500	Fax +61 2 9878 1055
Austria:	Vienna	Tel. +43 1 486 80 50 0	Fax +43 1 486 80 50 30
Belgium:	Groot Bijgaarden	Tel. +32 2 790 98 50	Fax +32 2 790 98 68
Canada:	Concord/Ontario	Tel. +1 800 248 0123	Fax +1 847 405 0164
Denmark:	Ballerup	Tel. +45 4454 0101	Fax +45 4454 0111
France:	Nanterre Cedex	Tel. +33 811 000 664	Fax +33 1 56 05 23 23
Germany:	Wetzlar	Tel. +49 64 41 29 40 00	Fax +49 64 41 29 41 55
Italy:	Milan	Tel. +39 02 574 861	Fax +39 02 574 03392
Japan:	Tokyo	Tel. +81 3 5421 2800	Fax +81 3 5421 2896
Korea:	Seoul	Tel. +82 2 514 65 43	Fax +82 2 514 65 48
Netherlands:	Rijswijk	Tel. +31 70 4132 100	Fax +31 70 4132 109
People's Rep. of China:	Hong Kong	Tel. +852 2564 6699	Fax +852 2564 4163
	Shanghai	Tel. +86 21 6387 6606	Fax +86 21 6387 6698
Portugal:	Lisbon	Tel. +351 21 388 9112	Fax +351 21 385 4668
Singapore		Tel. +65 6779 7823	Fax +65 6773 0628
Spain:	Barcelona	Tel. +34 93 494 95 30	Fax +34 93 494 95 32
Sweden:	Kista	Tel. +46 8 625 45 45	Fax +46 8 625 45 10
Switzerland:	Heerbrugg	Tel. +41 71 726 34 34	Fax +41 71 726 34 44
United Kingdom:	Milton Keynes	Tel. +44 800 298 2344	Fax +44 1908 246312
USA:	Buffalo Grove/Illinois	Tel. +1 800 248 0123	Fax +1 847 405 0164

and representatives in more than 100 countries

