



Leica TCS LSI

The World's First Super Zoom Confocal!
Technical Documentation

Specifications

Scan Head	Scanner	Method	true confocal	
		Confocal channels	1	
		Scanner	galvo (x,y)	
		Sequential scan	yes	
		Channel multiplexing	1 – 8 sequential	
	Resolution	Range (min – max)	[pixel]	128 ² – 2048 ²
		Scan formats	[pixel]	128, 256, 512, 1024, 2048
		Image depth	[bit]	8 or 12, switchable
	Spectral Detection	Spectral detection		yes
		Type		continuously variable
		Spectral resolution	[nm]	10 nm
		Bandwidth	[nm]	430 – 750
	Detector	Detector		1
		Detector type		ultra high dynamic PMT
		Detector connection		direct
	Confocal Zoom	Zoom type		continuously variable
Zoom range		[x]	1x – 16x	
Zoom increment			0.1	
Optical Zoom	Zoom type		motorized/manual	
	Zoom range, Z16 APO/Z6 APO	[x]	0.57 – 9.2x/0.57 – 3.6x	
	Zoom increment		continuously / variable	
Pinhole	Pinhole type		motorized, variable	
	Range (min – max)	[µm]	35 – 600	
	Pinhole adjustment	[%]	0 – 100	
	Control		automated via GUI	
Beam Splitter	Type		high performance dichroics	
	Beam splitter wavelength	[nm]	405/532; 488/635	
	ND-Splitter for TLD	[%]	RT 30/70	
	Beam splitter change		automated	
Scan Modes	2D: x-y, x-z, x-t		yes	
	3D: x-y-z, x-y-t, x-y-λ, x-z-y		yes	
	4D: x-y-z-t		yes	
Speed	Speed mode		uni/bi-directional	
	Line speed	[Hz]	800	
	Line speed range	[Hz]	400, 600, 800	
	Max@ 128 ²	[f/s]	6.0	
	standard@ 512 ²	[f/s]	2.0	
	Min@2048 ²	[f/s]	0.36	
Motorized Zooms	Z16 APO A	Type	Z16 APO A	
		Magnification range	0.57 – 9.2	
		NA-range depending on magnification	[nA] 0.017 – 0.112 (with Obj. 1x)	
		Zoom positioning	continuously/UMC controlled	
		Fine focus optics	yes motorized	
		Diaphragma	yes motorized	
	Z6 APO A	Type	Z6 APO A	
		Magnification range	0.57 – 3,6	
		NA-range depending on magnification	[nA] 0.02 – 0.117 (with Obj. 1x)	
		Positioning	continuously UMC controlled	
		fine focus optics	yes motorized	
		Diaphragma	yes motorized	

Supply Unit	Laser	Laser type		solid state	
		Number of lasers		max 4	
	Laser options	[nm]		405, 488, 532, 635	
	Excitation attenuation			AOTF	
	Excitation attenuation control			automated	
	Range	[%]		0 – 100	
	Typ.power@Zoom 1, Obj 1x	[mW]		> 1	
	Power depending on optical zoom			yes	
	Computer	PC integration			yes
		Processor			Intel Pentium M
		HD-size	[GB]		120
	Interface	Operating system			XP-embedded
		Ethernet			1
		USB			4
		Fire wire			2
		Parallel			1
	Mouse	Serial			1
		Type			optical
	Monitor	Monitor resolution	[Pixel]		1280x1024
		Monitor size			19"
Power Supply	Power supply integration			yes	
	Type			autoselect	
	Voltage range	[V]		100 – 240	
	Power consumption			800 VA	
Z-Drive	Micro	Z-focus		Galvanometer-stage	
		Z-stepsize	[@500 µm range]	10 nm	
		Max Z-Range		1500 µm	
Microscope	Types	Microscope upright		Base CSQ	
		Focussing drive		motorized	
		Max travel range	[mm]		150
	Laser Safety Chamber	Door angle, open	[°]		180
		Material			lasersafe acrylic glas
		Inside height from ground plate	[mm]		180
		Depth min	[mm]		420
		Width inside	[mm]		555
		Heating unit connection			yes
		Heating unit connection, diameter	[mm]		100
		Holder for micromanipulators			integrated
	Transmitted light	Intensity control			potentiometer
	Objectives	Macro-Objective	Magnification		1x, 2x, 5x
			NA	[mm]	variable
Working distance			[mm]	97/39/19	
Micro-Objectives ACS		Magnification			10x, 40x, 63x
		NA			0.3; 0.8; 1.15
		Working distance			0.3; 0.16; 0.15
Software	LAS AF LSI	Basic aquisition license		system included	
		Life Data Mode		optional	
		Colocalisation		optional	
		Deconvolution		optional	
		3D-visualization		optional	
		Dye Finder		optional	
Export formats			lif, tiff, avi, jpg		
Environment	Ratings	Humidity	[%]	10 – 80	
		Operating temperature	[°C]	18 – 30	
		Guaranteed stability	[°C/h]	23 +/- 2	


Zoom/Objective Combinations:

Note: Conditions recommended for confocal applications = green, acceptable = light green, not recommended = white.


Macro							Micro						
Objective	Z16 APO A			Z 6 APO A			Z16 APO A			Z 6 APO A			
Zoom Type	1x	2x	5x	1x	2x	5x	10x	40x	63x	10x	40x	63x	
Objective													
Free Working Distance [mm]	97	39	19	97	39	19	0.30	0.16	0.15	0.30	0.16	0.15	
Zoom													
Label	Objective Magnification Factor [x]						Objective Magnification Factor [x]						
0.57	0.7	1.4	3.5	0.7	1.4	3.5	0.57	3.5	14.1	22.3	3.5	14.1	22.3
0.8	1.0	2.0	5.0	1.0	2.0	5.0	0.8	5.0	19.8	31.2	5.0	19.8	31.2
1	1.2	2.5	6.2	1.2	2.5	6.2	1	6.2	24.8	39.1	6.2	24.8	39.1
1.25	1.6	3.1	7.8	1.6	3.1	7.8	1.25	7.8	31.0	48.8	7.8	31.0	48.8
1.6	2.0	4.0	9.9	2.0	4.0	9.9	1.6	9.9	39.7	62.5	9.9	39.7	62.5
2	2.5	5.0	12.4	2.5	5.0	12.4	2	12.4	49.6	78.1	12.4	49.6	78.1
2.5	3.1	6.2	15.5	3.1	6.2	15.5	2.5	15.5	62.0	97.7	15.5	62.0	97.7
3.2	4.0	7.9	19.8	4.0	7.9	19.8	3.2	19.8	79.4	125.0	19.8	79.4	125.0
3.6	4.5	8.9	22.3	4.5	8.9	22.3	3.6	22.3	89.3	140.6	22.3	89.3	140.6
4.6	5.7	11.4	28.5				4.6	28.5	114.1	179.7			
5	6.2	12.4	31.0				5	31.0	124.0	195.3			
6.3	7.8	15.6	39.1				6.3	39.1	156.2	246.1			
8	9.9	19.8	49.6				8	49.6	198.4	312.5			
9.23	11.4	22.9	57.2				9.23	57.2	228.9	360.5			

Field of View [mm]							Field of View [mm]						
0.57	22.0	11.0	4.4	22.0	11.0	4.4	0.57	4.40	1.10	0.70	4.40	1.10	0.70
0.8	15.7	7.8	3.1	15.7	7.84	3.14	0.8	3.14	0.78	0.50	3.14	0.78	0.50
1	12.5	6.3	2.5	12.5	6.27	2.51	1	2.51	0.63	0.40	2.51	0.63	0.40
1.25	10.0	5.0	2.0	10.0	5.02	2.01	1.25	2.01	0.50	0.32	2.01	0.50	0.32
1.6	7.8	3.9	1.6	7.8	3.92	1.57	1.6	1.57	0.39	0.25	1.57	0.39	0.25
2	6.3	3.1	1.3	6.3	3.14	1.25	2	1.25	0.31	0.20	1.25	0.31	0.20
2.5	5.0	2.5	1.0	5.0	2.51	1.00	2.5	1.00	0.25	0.16	1.00	0.25	0.16
3.2	3.9	2.0	0.8	3.9	1.96	0.78	3.2	0.78	0.20	0.12	0.78	0.20	0.12
3.6	3.5	1.7	0.7	3.5	1.74	0.70	3.6	0.70	0.17	0.11	0.70	0.17	0.11
4.6	2.7	1.4	0.5				4.6	0.55	0.14	0.09			
5	2.5	1.3	0.5				5	0.50	0.13	0.08			
6.3	2.0	1.0	0.4				6.3	0.40	0.10	0.06			
8	1.6	0.8	0.3				8	0.31	0.08	0.05			
9.23	1.4	0.7	0.3				9.23	0.27	0.07	0.04			

Please note: an increase of NA and resolution is achieved up to zoom value 6.3.









- 1 Confocal scanhead
- 2 Supply unit
- 3 Laser safety chamber
- 4 Wing doors
- 5 Z-zoom, motorized
- 6 SuperZ Galvo stage
- 7 Motor focus drive
- 8 DFC camera option
- 9 Anti-vibration table, passive
- 10 Monitors
- 11 Mouse
- 12 Keyboard
- 13 UMC control
- 14 XY-stage control
- 15 Heat pipe adapter
- 16 EL6000 fluorescence illumination
- 17 Computer table



System Overview Leica TCS LSI

Solid State Laser:




-  R: SS 10 mW 635 nm
-  G: SS 10 mW 532 nm
-  B: SS 10 mW 488 nm
-  V: SS 10 mW 405 nm



15 6901 101 Leica TCS LSI with Base CSQ VIS

Additional options:

Digital Camera Kit:

-  **12 7300 42**
Leica DFC300 FX
-  **12 7300 43**
Leica DFC350 FX¹
-  **15 6905 300**
LSI Camera Kit for DFC300/350, Including: Software, C-Mount Adapter 1x, Cable

1) 04/2008: DFC 360 FX

SuperZ Galvo Stage Inserts:





-  **15 6604 410**
Stage Insert for SuperZ, rotatable
-  **15 5935 411**
Stage Insert for SuperZ, universal
-  **15 6604 412**
Stage Insert for SuperZ, basic (system included)

Plate Inserts:

- 15 6604 420** Plate Insert - Microtiter Plate
- 15 6604 421** Plate Insert - Willco Dish 35 mm
- 15 6604 422** Plate Insert - Willco Dish 50 mm
- 15 6604 423** Plate Insert - MatTek Dish 35 mm
- 15 6604 424** Plate Insert - Petri Dish 36 mm
- 15 6604 425** Plate Insert - Petri Dish 39 mm

Climate Control:

- 15 6905 400** Leica LSI Climate Kit
Including Tubing, Sensor Heating Control Unit
- 15 6905 600** TCS LSI CO₂-Control Set 110V/230V
Including: SuperZ-CO₂-Cover
CO₂-Controller, Humidifier
- 15 6905 602** Universal Gas Cover for SuperZ with In/Outlet



15 6901 110 Leica TCS LSI with Base CSQ V-VIS

Zoom-Optics:

- ① **15 6904 602**
Zoom 6x APO motorized Kit
- ② **15 6904 603**
Zoom 16x APO motorized Kit

Macro-Objectives for Zoom Optics

- ② **15 6904 610**
Plan APO 1x Objective
- 15 6904 611**
Plan APO 2x Objective
- 15 6904 612**
Plan APO 5x Objective
- 15 6904 620** ⑮
Objective Zoom Adapter

Micro-Objectives for Z16 APO (A)

- ⑮ **15 6904 630**
10x Objective
- 15 6904 631**
40x Objective
- 15 6904 632**
63x W Objective
- 156904633**
63x Oil Objective

Manual XY-Stage Kits:

- 15 6905 210**
ISO pro Stage manual Kit
incl. Cover and Z-adapter
- 15 6905 211**
LSI 3-Plate manual Stage
Kit incl. Cover

LSI:

- 15 6906 100**
Microscope Table LSI,
passive
- 156906102**
Microscope Table LSI,
active
- 15 6906 101**
Computer Table

Motorized XY-Stage Kits:

- ④ **15 6905 200**
ISO pro Stage mot. Kit
with Drive Control, Joy
Stick and Z-adapter
- 15 6905 201**
LSI 3-Plate mot. Kit
Including
CTR 6500 Drive Control,
Cover Plate
- 15 6905 202**
LSI mot. Scanning
Stage Kit²
Including
CTR 6500 Drive Control
Base Adapters

Software LAS AF - TCS LSI

- ⑪ **15 69 02 202**
LAS AF LSI Life Data Mode
- 15 6902 203**
LAS AF LSI Co-Localization
- 15 6902 204**
LAS AF LSI Deconvolution
- 15 6902 205**
LAS AF LSI 3D Visualization
- 15 6902 208**
LAS AF LSI Dye Finder
- 15 6902 211**
LAS AF LSI Camera Advanced


Monitor-Option

- ⑮ **15 6905 700**
Second 19"-TFT Monitor

Room Requirements



visible radiation:



Power Supply:	Power supply integration	yes
	Type	autoselect
	Voltage range	[V] 100 – 240
	Power consumption	[VA] 800
	Independent circuits	[no.] 1
	Frequency	[Hz] 50/60
	Fuse: standard	[A] 10

Note: The optimal optical performance can only be achieved on stable room floors. Concrete floors are required. Other like e.g. wooden floors are not suitable.

Environment	Humidity	[%]	10 – 80
	Operating temperature	[°C]	18 – 30
	Guaranteed stability		23°C +/- 2°C

Load Capacity and Weight	Confocal unit, max.	[kg]	75
	Microscope, max.	[kg]	45
	System	[kg]	90
	Static floor load	[kg/m ²]	200

