




1. Data Transfer Hub



2. USB



3. Data Transfer Hub



4.



60 sec.



5. App Store




6. WiFi





7. WiFi



8. Leica DMshare



9. Photo



10. Setup



WiFi Access Point

Camera
Live image size:

Admin
 Require PIN

Set a PIN to restrict access to this page.

Networking
 Built-in WiFi Access Point
Password:
 Join existing WiFi

Local WiFi

Require PIN

Set a PIN to restrict access to this page.

Networking
 Built-in WiFi Access Point
 Join existing WiFi
Name (SSID):
Password:

AP: The HUB can act as a WiFi Access Point allowing the iPad connect directly - no infrastructure required. Though you will need to connect your iPad to the HUBs own WiFi network (DMshare-PT0001) which will prevent internet access.

Existing Network: you can connect the DMshare Hub to your existing WiFi network - the one your iPad is already connected to - so you don't need to switch back and forth for internet access.



Leica DMshare 2.0

Quick Start Guide