



ARTICLE DESCRIPTIONS



Cables
 11 600 269, FW-b-b
 11 600 254, FW-a-b
 12 730 211, HDMI-cable
 12 730 530, USB3-cable

PC Cards
 12 730 446, FW-b, notebook kit
 12 730 447, FW-b, PCI express
 12 730 495, USB3, PCI express

Leica Software PC



Name	MPixel	Sensor	Size	Pixel-Size	Cooling	Interface	C-Mount*
MC170 HD	5	CMOS-c	1/2.3"	2.35 µm	No	USB2+HD	0.4x-0.5x
MC190 HD	10	CMOS-c	1/2.3"	1.67 µm	No	USB2+HD	0.4x-0.5x
DMC2900	3	CMOS-c	1/2"	3.2 µm	No	USB3	0.5x
DMC4500	5	CCD-c	2/3"	3.4 µm	No	USB3	0.63x
DFC450 C	5	CCD-c	2/3"	3.4 µm	Yes	FWb + Sync	0.5x
DFC7000 T	2.8	CCD-c	2/3"	4.54 µm	Yes	USB3	0.63x
DFC9000 GT	4.2	sCMOS	19 mm	6.5 µm	Yes	USB3 + Camlink	1x
DFC9000 GTC	4.2	sCMOS	19 mm	6.5 µm	Yes	USB3 + Camlink	1x
DMC5400	20	CMOS-c	1"	2.4 µm	No	USB3	1x
DMC6200	2.3/20	CMOS-c	1.2"	5.86 µm	No	USB3	1x
DFC3000 G	1.3	CCD-m	1/3"	3.75 µm	No	USB3 + Sync	0.32x-0.4x

-c: color / -m: monochrome / * recommended


Compound 2017/V.1

Digital camera systems
 12 730 465 Leica MC170 HD camera kit
 12 730 520 Leica MC190 HD camera kit
 12 730 517 Leica DMC4500 camera kit
 12 730 466 Leica DMC2900 camera kit
 12 730 412 Leica DFC450 C camera kit
 12 730 531 Leica DMC5400 camera kit
 12 730 532 Leica DMC6200 camera kit

11 547 005 Leica DFC3000 G camera kit
 11 547 106 Leica DFC7000T camera kit
 11 547 110 Leica DFC9000 GT camera kit
 11 547 007 Leica DFC9000 GTC camera kit


Digital camera system accessories
 11 600 269 FireWire cable, FW-b-b, 2.5 m, 9-pin to 9-pin
 12 730 211 HDMI-cable, 3 m, HDMI 2x standard-A
 12 730 530 USB3-cable, 3 m standard-A to micro-B
 12 730 447 FW-b, PCI express card
 12 730 495 USB3, PCI express card

Phototubes and camera adapter
 10 447 319 Adapter T2, Canon EOS
 10 446 175 2.5x DSLR tube with T2 thread (24x36 mm sensor)
 11 621 500 MPV adapter, HC/delta
 11 541 038 Vario adapter for C-mount 0.55x to 1.1x – Delta
 11 541 510 C-mount adapter HC 1x, for 1" digital cameras
 11 541 543 C-mount adapter HC 0.7x, for 2/3" digital cameras
 11 541 544 C-mount adapter HC 0.55x, for 1/2" digital cameras
 10 450 622 C-mount adapter HC 0.4x, for 1/2" digital cameras
 11 541 512 C-mount adapter HC 0.35x, for 1/3" digital cameras




Connection to main diagram next page

ARTICLE DESCRIPTIONS



Cables
 11 600 269, FW-b-b
 11 600 254, FW-a-b
 12 730 211, HDMI-cable
 12 730 530, USB3-cable

PC Cards
 12 730 446, FW-b, notebook kit
 12 730 447, FW-b, PCI express
 12 730 495, USB3, PCI express




MC170 HD
 12 730 465
MC190 HD
 12 730 520

DMC2900
 12 730 466
DMC4500
 12 730 517

DFC3000 G
 11 547 005
DFC450 C
 12 730 412

DMC5400
 12 730 531
DMC6200
 12 730 532

DFC7000 T
 11 547 106
DFC9000 GT
 11 547 110
DFC9000 GTC
 11 547 007



10 445 928
0.32x

10 450 528
0.5x

10 447 367
0.63x

10 446 307
0.8x

10 450 829
1x

HC
 11 541 510 1.0x
 11 541 543 0.7x
 11 541 544 0.55x
 10 450 622 0.4x


Delta
 10 446 175 2.5x
 11 541 006 1.0x
 11 541 007 0.63x
 11 541 016 0.5x



10 450 317
1.0x



10 445 930
1.0x



10 447 436
1.6x

1
**Connection to main diagram
 next page**

Name	MPixel	Sensor	Size	Pixel-Size	Cooling	Interface	C-Mount*
MC170 HD	5	CMOS-c	1/2.3"	2.35 µm	No	USB2+HD	0.4x-0.5x
MC190 HD	10	CMOS-c	1/2.3"	1.67 µm	No	USB2+HD	0.4x-0.5x
DMC2900	3	CMOS-c	1/2"	3.2 µm	No	USB3	0.5x
DMC4500	5	CCD-c	2/3"	3.4 µm	No	USB3	0.63x
DFC450 C	5	CCD-c	2/3"	3.4 µm	Yes	FWb+ Sync	0.5x
DFC7000 T	2.8	CCD-c	2/3"	4.54 µm	Yes	USB3	0.63x
DFC9000 GT	4.2	sCMOS	19 mm	6.5 µm	Yes	USB3 + Camlink	1x
DFC9000 GTC	4.2	sCMOS	19 mm	6.5 µm	Yes	USB3 + Camlink	1x
DMC5400	20	CMOS-c	1"	2.4 µm	No	USB3	1x
DMC6200	2.3/ 20	CMOS-c	1.2"	5.86 µm	No	USB3	1x
DFC3000 G	1.3	CCD-m	1/3"	3.75 µm	No	USB3+ Sync	0.32x-0.4x

Digital camera systems

- 12 730 465 Leica MC170 HD camera kit
- 12 730 520 Leica MC190 HD camera kit
- 12 730 517 Leica DMC4500 camera kit
- 12 730 466 Leica DMC2900 camera kit
- 12 730 412 Leica DFC450 C camera kit
- 12 730 531 Leica DMC5400 Kamera Kit
- 12 730 532 Leica DMC6200 Kamera-Kit

- 11 547 005 Leica DFC3000 G camera kit
- 11 547 106 Leica DFC7000T camera kit
- 11 547 110 Leica DFC9000 GT camera kit
- 11 547 007 Leica DFC9000 GTC camera kit

- 12 730 521 Leica IC90 E camera (incl. USB-cable and Leica software) not recommended for MZ10 F / M165 FC / M205 FA
- 12 730 528 Leica stand alone kit (HDMI cable, SD card, USB power supply, remote control RC3, ethernet cable) for IC90 E
- 12 730 229 Hand or foot switch with 2 m cable for IC90 E

Digitale camera system accessories

- 11 600 269 FireWire cable, FW-b-b, 2.5 m, 9-pin to 9-pin
- 12 730 211 HDMI-cable, 3 m, HDMI 2x standard-A
- 12 730 530 USB3-cable, 3 m standard-A to micro-B
- 12 730 447 FW-b, PCI express card
- 12 730 495 USB3, PCI express card

Phototubes and C-mounts

- 10 447 319 Adapter T2, Canon EOS
- 10 447 436 1.6x DSLR tube with T2 thread (APS-C sensor)
- 10 446 175 2.5x DSLR tube with T2 thread (24x36 mm sensor)
- 10 445 928 0.32x video objektive with C-mount for 1/3" digital cameras
- 10 450 528 0.5x video objektive with C-mount for 1/2" digital cameras
- 10 447 367 0.63x video objektive with C-mount for 2/3" digital cameras
- 10 446 307 0.8x video objektive with C-mount for 2/3" digital cameras
- 10 445 930 1.0x video-/photo objektive
- 10 450 317 1.0x C-mount adapter for HC