

From Eye to Insight



Leica Microsystems (Schweiz) AG | Max Schmidheiny-Strasse 201 | 9435 Heerbrugg (Switzerland)

To whom it may concern,

19th August 2024


Antibacterial Activity Evaluation Test declaration

Proof of antibacterial evaluation test according to JIS Z 2801:2000

Coating of below mentioned products pass antibacterial activity evaluation test according to JIS Z 2801:2000 with a value of 3.7 for the antibacterial activity (Which corresponds to a killing ratio of >99.97 %) for Staphylococcus aureus (ATCC 6538P) The test rapport confirms an effectivity of >3.5 which means >99.97% killing ratio.

Order Nr.	Product
13 613 215	DM500 RH
13 613 216	DM500 LH
13 613 017	DM 750/5K LH
13 613 015	DM 750/4K LH
13 613 016	DM 750/5 LH
13 613 014	DM 750/4 LH
13 613 013	DM 750/5K
13 613 011	DM 750/4K
13 613 012	DM 750/5
13 613 010	DM 750/4
13 613 615	DM 750 P

The antibacterial activity evaluation test according to JIS Z 2801:2000 was done by Company Sanitized AG, test report 2008-0881.

DocuSigned by:

12B020BA97F540E...

Dr. Jochen Langer
Senior Manager Global Product Compliance