



Product Service

# CERTIFICATE


No. Z1 17 05 56257 161

**Holder of Certificate:** Leica Microsystems CMS GmbHErnst-Leitz-Straße 17-37  
35578 Wetzlar  
GERMANY**Certification Mark:****Product:** Laboratory equipment  
(microscope system)

The product was tested on a voluntary basis and complies with the essential requirements. The certification mark shown above can be affixed on the product. It is not permitted to alter the certification mark in any way. In addition the certification holder must not transfer the certificate to third parties. See also notes overleaf.

**Test report no.:** 028-713099377-000**Valid until:** 2022-05-25**Date,** 2017-05-30  
( Ralph Fischer )

Page 1 of 2





Product Service

**CERTIFICATE**  
**No. Z1 17 05 56257 161**

**Model(s):** **The microscope system consists of:  
WSU, Pulsed Laser Unit 355 nm,  
Infinity scanner, Supply Unit Scanhead**

**Brand Name:** **Leica**

**Parameters:**

WSU	
Rated voltage:	100-240 V
Rated frequency:	50/60 Hz
Rated power:	250 VA
Protection class:	I
Supply Unit Scanhead (supplied Infinity Scanner)	
Rated voltage:	100-240 V
Rated frequency:	50/60 Hz
Rated power:	150 VA
Protection class:	I
Pulsed Laser Unit 355 nm (Ablation laser)	
Rated voltage:	12 VDC (supplied via external limited power supply unit)
Protection class:	III
Laser class:	1
Protection class:	I

**Tested according to:** EN 61010-1:2010  
EN 60825-1:2014

**Production Facility(ies):** 56257

Page 2 of 2