



America

CERTIFICATE

No. U14 056257 0310 Rev. 00

Holder of Certificate: **Leica Microsystems CMS GmbH**
Ernst-Leitz-Straße 17-37
35578 Wetzlar
GERMANY

Certification Mark:



C US

Product: **Laboratory equipment
(Ultramicrotome)**

Brand Name: **Leica**

Model(s): UC ENUITY system consist of:
UC ENUITY Cryo-Chamber, UC ENUITY Power unit,
UC ENUITY LN2 Pump

Parameters:

Rated voltage:	100-240 VAC
Rated frequency:	50/60 Hz
Rated power:	350 VA
Protection class:	I

Tested according to:

- UL 61010-1:2012/R:2019-07
- UL 61010-2-010:2019
- UL 61010-2-011:2021
- CSA C22.2 No. 61010-1:2012/A1:2018-11
- CSA C22.2 No. 61010-2-010:2019
- CSA C22.2 No. 61010-2-011:2019

The product was tested on a voluntary basis and complies with the USA and Canadian standards stated on this certificate. It can be marked with the certification mark above. The certification mark shown above can be affixed on the product. It is not permitted to alter the certification mark in any way. In addition, the certification holder must not transfer the certificate to third parties.

Test report no.: 713307996-01

Date, 2023-11-25

Siemon

(Thorsten Siemon)